

GENEzol™ 96 Well TriRNA Pure Kit

96GZX02 (2 x 96 well plates/kit)

96GZX04 (4 x 96 well plates/kit)

96GZX10 (10 x 96 well plates/kit)

Advantages

Sample: up to 2×10^9 cultured cells, 10-20 mg of animal tissue, up to 150 μ l of body fluids, up to 1×10^9 bacteria cells, 10-20 mg of plant tissue per well

Binding Capacity: 50 μ g of RNA per well

Format: GENEzol™ Reagent combined with Presto™ RNA 96 Well Binding Plate

Operation Time: 60 minutes

Elution Volume: 60-80 μ l (dead volume: 20-25 μ l)

Kit Storage: dry at room temperature (15-25°C)

Table of Contents

| | |
|--------------------------------------|----|
| Introduction..... | 2 |
| Quality Control..... | 2 |
| Kit Components..... | 2 |
| Safety Measures..... | 2 |
| Quick Protocol Diagram..... | 3 |
| Sample Homogenization and Lysis..... | 4 |
| Centrifuge Protocol Procedure..... | 5 |
| Vacuum Protocol Procedure..... | 7 |
| Troubleshooting..... | 9 |
| Test Data..... | 10 |
| Related Products..... | 11 |

Introduction

The GENEzol™ 96 Well TriRNA Pure Kit is a phenol and guanidine isothiocyanate plus 96 well RNA binding plate system for high-throughput purification of high-quality total RNA from a variety of samples. Initially, samples are homogenized in GENEzol™ Reagent without chloroform phase separation or isopropanol RNA precipitation. Following sample homogenization, simply bind, wash and elute the high-quality, total RNA in RNase-free Water and use in a variety of sensitive downstream applications.

Quality Control

The quality of the GENEzol™ 96 Well TriRNA Pure Kit is tested on a lot-to-lot basis according to Geneaid's ISO-certified quality management system. 10 µl from a 50 µl eluate of purified RNA is analyzed by electrophoresis on a 0.8% agarose gel.

Kit Components

| Component | 96GZX02 | 96GZX04 | 96GZX10 |
|---|---------------------------|---------------------------|-----------------------------|
| GENEzol™ Reagent | 100 ml | 200 ml | 250 ml x 2 |
| Pre-Wash Buffer ¹ (Add Ethanol) | 70 ml (30 ml) | 70 ml x 2 (30 ml x 2) | 175 ml x 2 (75 ml x 2) |
| Wash Buffer ² (Add Ethanol) | 50 ml x 2 (200 ml x 2) | 50 ml x 4 (200 ml x 4) | 50 ml x 10 (200 ml x 10) |
| RNase-Free Water | 30 ml | 15 ml x 1 30 ml x 1 | 30 ml x 3 |
| Presto™ RNA 96 Well Binding Plates | 2 | 4 | 10 |
| Microtubes (Racked) | 2 | 2 | 2 |
| Microtubes (8-strip) | 12 x 2 | 12 x 6 | 12 x 18 |
| Caps for Microtubes (8-strip) | 48 | 48 x 2 | 48 x 5 |
| 96 Deep Well Plates ³ | 2 | 2 | 2 |
| 0.35 ml Collection Plates | 2 | 4 | 10 |
| Adhesive Film | 2 | 4 | 10 |


^{1,2}Add absolute ethanol (see the bottle label for volume) to Pre-Wash Buffer and Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

³96 Deep Well Plates are reusable. After use, rinse the plate with water then incubate in 0.4M HCl for 1 minute at room temperature. Wash the plate thoroughly with ddH₂O. The plate can be autoclaved after being washed.



GENEzol™ Reagent contains phenol and guanidine isothiocyanate. During operation, always work in a fume hood, always wear a lab coat, disposable gloves, protective goggles and (anti-fog) procedure mask. Disposable/non-disposable glassware, plasticware and automatic pipettes should be sterile (RNase-free) and used only for RNA procedures.

Quick Protocol Diagram



GENEzol™ 96 Well TriRNA Kit Protocol

Please read the entire instruction manual prior to starting the Protocol Procedure.

IMPORTANT BEFORE USE!

1. Add absolute ethanol (see the bottle label for volume) to Pre-Wash Buffer and Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

2. 96 Deep Well Plates are reusable. After use, rinse the plate with water then incubate in 0.4M HCl for 1 minute at room temperature. Wash the plate thoroughly with ddH₂O. The plate can be autoclaved after being washed.

Additional Requirements

Centrifuge with microplate buckets/vacuum manifold, additional 96 Deep Well Plate, absolute ethanol, steal or ceramic beads, TissueLyser or mixer mill for animal/plant tissue. Lysozyme (LY420) and Bacteria Lysis Buffer (BLB00030) for bacteria, (optional) RNase-Free DNase I Set (DNS50/100/300) can be purchased directly from Geneaid.

Sample Homogenization and Lysis (Centrifuge and Vacuum)

| Sample | Procedure |
|---------------------------|---|
| Adherent Cultured Cells | <ol style="list-style-type: none"> 1. Remove the culture medium from the culture dish. 2. Directly add 500 µl of GENEzol™ Reagent to the culture dish (up to 2 x 10⁶ cells) then lyse the cells directly in the culture dish by pipetting several times. 3. Incubate for 5 minutes at room temperature. 4. Transfer the sample to each tube of Microtubes (Racked), being careful not to touch the rims with pipette tips. |
| Suspension Cultured Cells | <ol style="list-style-type: none"> 1. Transfer cells (up to 2 x 10⁶) to each tube of Microtubes (Racked), being careful not to touch the rims with pipette tips. 2. Seal the tubes with Caps for Microtubes then cover rack with the plastic cover. 3. Centrifuge at 300 x g for 5 minutes. 4. Remove caps then completely remove the culture medium. 5. Add 500 µl of GENEzol™ Reagent into each tube then re-suspend the cell pellet by pipette. 6. Incubate for 5 minutes at room temperature. |
| Animal and Plant Tissue | <ol style="list-style-type: none"> 1. Add steal beads or ceramic beads (RNase-free) into each tube of Microtubes (Racked). 2. Transfer 10-20 mg of animal or plant tissue and 500 µl of GENEzol™ Reagent into each tube, being careful not to touch the rims with pipette tips. 3. Seal the tubes with Caps for Microtubes then cover rack with the plastic cover. 4. Homogenize samples using a TissueLyser or mixer mill then incubate for 5 minutes at room temperature. 5. Centrifuge at 5,000-6,000 x g for 5 minutes to remove cell debris then transfer the clear supernatant to each tube of a new Microtubes (Racked). |

| | |
|---|---|
| Body Fluids (blood, buffy coat, plasma, serum) | <ol style="list-style-type: none"> 1. Transfer up to 150 µl of liquid sample to each tube of Microtubes (Racked) then add 3 volumes of GENEzol™ Reagent per 1 volume of sample (3:1) and mix well by pipette, being careful not to touch the rims with pipette tips. 2. Seal the tubes with Caps for Microtubes then cover rack with the plastic cover. 3. Incubate for 5 minutes at room temperature. 4. Centrifuge at 5,000-6,000 x g for 5 minutes to remove cell debris then transfer the clear supernatant to each tube of new Microtubes (Racked). |
| Bacteria | <ol style="list-style-type: none"> 1. Transfer bacterial cells (up to 1 x 10⁹) to each tube of Microtubes (Racked), being careful not to touch the rims with pipette tips. 2. Seal the tubes with Caps for Microtubes then cover rack with the plastic cover. 3. Centrifuge at 5,000-6,000 x g for 3 minutes. 4. Remove caps then completely remove the culture medium. 5. Transfer 100 mg of Lysozyme powder to a 15 ml centrifuge tube (RNase-free) containing 10 ml of Bacteria Lysis Buffer then vortex until the Lysozyme powder is completely dissolved. 6. Add 100 µl of Bacteria Lysis Buffer containing Lysozyme into each tube then re-suspend the cell pellet by pipette. 7. Incubate for 5 minutes at room temperature. 8. Add 500 µl of GENEzol™ Reagent into each tube then mix well by pipette. 9. Incubate for 5 minutes at room temperature. |

Centrifuge Protocol Procedure

2. RNA Binding

Add 1 volume of absolute ethanol directly to 1 volume of sample mixture (1:1) into each tube of **Microtubes (Racked)**, being careful not to touch the rims with pipette tips. Seal the tubes with new caps, and cover the rack with the plastic cover. Mix by shaking vigorously for 15-30 seconds. Centrifuge at 1,000 x g for 30 seconds to collect any lysate from the caps. Place the **Presto™ 96 Well RNA Binding Plate on a 96 Deep Well Plate**. Remove caps then transfer **500 µl of sample lysate** into each well of the **Presto™ 96 Well RNA Binding Plate**, being careful not to get any lysate on the the rims of the wells. Centrifuge the **Presto™ 96 Well RNA Binding Plate** and **96 Deep Well Plate** together at 3,000 x g for 5 minutes. Transfer the remaining sample lysate into each well of **Presto™ 96 Well RNA Binding Plate**, being careful not to get any lysate on the the rims of the wells. Centrifuge the **Presto™ 96 Well RNA Binding Plate** and **96 Deep Well Plate** together at 3,000 x g for 5 minutes. Discard the flow-through then place the **Presto™ 96 Well RNA Binding Plate** back on the **96 Deep Well Plate**.

Optional In-Column DNase I Treatment

Add 400 μ l of Wash Buffer (make sure ethanol was added) into each well of the Presto™ 96 Well RNA Binding Plate, being careful not to get any buffer on the the rims of the wells. Centrifuge the Presto™ 96 Well RNA Binding Plate and 96 Deep Well Plate together at 3,000 x g for 5 minutes. Discard the flow-through. Place the Presto™ 96 Well RNA Binding Plate back on the 96 Deep Well Plate.

Prepare DNase I solution in a 15 ml centrifuge tube (RNase-free) as follows:

| | |
|-------------------------|--------|
| DNase I (2U/ μ l) | 0.5 ml |
| DNase I Reaction Buffer | 4.5 ml |
| Total Volume | 5 ml |

Gently pipette the DNase I solution to mix (DO NOT vortex) then add 50 μ l of DNase I solution into the CENTER of each well of the Presto™ 96 Well RNA Binding Plate. Incubate for 15 minutes at room temperature (20-30°C) then proceed with RNA Wash.

3. Wash

Add 400 μ l of Pre-Wash Buffer (make sure ethanol was added) to each well of the Presto™ 96 Well RNA Binding Plate then centrifuge the Presto™ 96 Well RNA Binding Plate and 96 Deep Well Plate together at 3,000 x g for 5 minutes. Discard the flow-through then place the Presto™ 96 Well RNA Binding Plate back on the 96 Deep Well Plate. Add 700 μ l of Wash Buffer (make sure ethanol was added) to each well of the Presto™ 96 Well RNA Binding Plate then centrifuge the Presto™ 96 Well RNA Binding Plate and 96 Deep Well Plate together at 3,000 x g for 5 minutes. Discard the flow-through then place the Presto™ 96 Well RNA Binding Plate back on the 96 Deep Well Plate. Wash the Presto™ 96 Well RNA Binding Plate again with 700 μ l of Wash Buffer. Discard the flow-through then place the Presto™ 96 Well RNA Binding Plate back on the 96 Deep Well Plate. Centrifuge the Presto™ 96 Well RNA Binding Plate and 96 Deep Well Plate together at 3,000 x g for 10 minutes to dry the membrane.

4. Elution

Remove the **Presto™ 96 Well RNA Binding Plate** from the **96 Deep Well Plate** then blot the nozzles on a clean absorbent paper towel to remove residual ethanol. Place the **Presto™ 96 Well RNA Binding Plate** on top of a **0.35 ml Collection Plate**. **Add 60-80 μ l of RNase-free Water** into the CENTER of each well of the **Presto™ 96 Well RNA Binding Plate**. Let stand for at least 2 minutes to ensure the RNase-free water is absorbed by the membrane. Centrifuge the **Presto™ 96 Well RNA Binding Plate** and **0.35 ml Collection Plate** together at 3,000 x g for 5 minutes. Seal the **0.35 ml Collection Plate** with **Adhesive Film** and store the purified RNA at -70°C.

Vacuum Protocol Procedure

2. Vacuum Manifold Preparation

Place the **96 Deep Well Plate** on the manifold base. Place the binding top plate on the manifold base. Place the **Presto™ 96 Well RNA Binding Plate** on the binding top plate aperture. Attach the vacuum manifold to a vacuum source.

3. RNA Binding

Add 1 volume of absolute ethanol directly to 1 volume of sample mixture (1:1) into each tube of **Microtubes (Racked)**, being careful not to touch the rims with pipette tips. Seal the tubes with new caps, and cover the rack with the plastic cover. Mix by shaking vigorously for 15-30 seconds. Centrifuge at 1,000 x g for 30 seconds to collect any lysate from the caps. Remove caps then transfer sample lysate (approx. 1 ml) into each well of the **Presto™ 96 Well RNA Binding Plate**, being careful not to get any lysate on the the rims of the wells. Seal unused wells of the **Presto™ 96 Well RNA Binding Plate** with **Adhesive Film**. Apply vacuum at 15 inches Hg until sample passes through the **Presto™ 96 Well RNA Binding Plate** then turn off the vacuum. Discard the flow-through and re-assemble manifold.

Optional In-Column DNase I Treatment

Add 400 µl of Wash Buffer (make sure ethanol was added) into each well of the Presto™ 96 Well RNA Binding Plate, being careful not to get any buffer on the the rims of the wells. Apply vacuum at 15 inches Hg until sample passes through the Presto™ 96 Well RNA Binding Plate then turn off the vacuum.

Prepare DNase I solution in a 15 ml centrifuge tube (RNase-free) as follows:

| | |
|-------------------------|--------|
| DNase I (2U/µl) | 0.5 ml |
| DNase I Reaction Buffer | 4.5 ml |
| Total Volume | 5 ml |

Gently pipette the DNase I solution to mix (DO NOT vortex) then add 50 µl of DNase I solution into the CENTER of each well of the Presto™ 96 Well RNA Binding Plate. Incubate for 15 minutes at room temperature (20-30°C) then proceed with RNA Wash.

4. Wash

Add 400 µl of Pre-Wash Buffer (make sure ethanol was added) to each well of the **Presto™ 96 Well RNA Binding Plate**. Apply vacuum at 15 inches Hg until **Pre-Wash Buffer** passes through the **Presto™ 96 Well RNA Binding Plate** (approx. 10 seconds). Turn off the vacuum. Discard the flow-through and re-assemble the manifold. **Add 1 ml of Wash Buffer (make sure ethanol was added)** to each well of the **Presto™ 96 Well RNA Binding Plate**. Apply vacuum at 15 inches Hg until **Wash Buffer** passes through the **Presto™ 96 Well RNA Binding Plate** (approx. 10 seconds). Turn off the vacuum. Discard the flow-through and re-assemble the manifold. Wash the **Presto™ 96 Well RNA Binding Plate** again with **1 ml of Wash Buffer**. Apply vacuum for an additional 10 minutes to dry the membrane then turn off the vacuum.

5. Elution

Remove the **Presto™ 96 Well RNA Binding Plate** from the binding top plate aperture then blot the nozzles on a clean absorbent paper towel to remove residual ethanol. Remove the **96 Deep Well Plate** from the manifold base then place the collection plate spacer on the manifold base. Place a **0.35 ml collection Plate** on top of the collection plate spacer. Place the binding top plate back on the manifold base then place the **Presto™ 96 Well RNA Binding Plate** back on the binding top plate aperture. **Add 60-80 µl of RNase-free Water** into the CENTER of each well of the **Presto™ 96 Well RNA Binding Plate**. Let stand for at least 2 minutes to ensure the **RNase-free Water** is absorbed by the membrane. Apply vacuum at 15 inches Hg for 5 minutes then turn off the vacuum. Seal the **0.35 ml Collection Plate** with **Adhesive Film** and store the purified RNA at -70°C.

Troubleshooting



Low Yield

Sample lysis or homogenization was incomplete.

Starting material should be reduced and completely dissolved in GENEzol™ Reagent.

Incorrect RNA elution.

Make sure RNase-free Water is added to the center of each well of the 96 Well RNA Binding Plate and is absorbed completely.

Precipitates may form during the RNA binding step after adding 1 volume of absolute ethanol to the sample mixture in GENEzol™ Reagent if too much sample material is used.

Reduce sample to half of the original amount.

Degraded RNA

Incorrect sample preparation and/or storage.

Process or freeze samples immediately after collection.

Incorrect storage temperature.

Extracted RNA should be stored at -70°C.

Low RNA A260/A280

Volume of GENEzol™ Reagent was insufficient for proper sample homogenization.

Volume of GENEzol™ Reagent is sample dependent and should be added according to the sample homogenization specifications.

Incomplete wash step.


Wash the 96 Well RNA Binding Plate with appropriate volume of ethanol added Wash Buffer.

Eluted RNA does not perform well in downstream applications

Residual ethanol contamination.

Following the wash step, dry the 96 Well RNA Binding Plate with additional centrifugation at 3,000 x g or with additional vacuum for 10 minutes to ensure the membrane is completely dry.

GENEzol™ 96 Well TriRNA Pure Kit Functional Test Data (HeLa Cells)



1 2 M 3 4

Figure 1. RNA was purified using the GENEzol™ 96 Well TriRNA Pure Kit in parallel to the similar product from competitor Z. 5×10^5 HeLa cells were homogenized using GENEzol™ Reagent and competitor Z tri reagent. RNA was then purified using the corresponding kits binding plate procedure. 10 μ l from a 50 μ l eluate of purified RNA was analyzed by electrophoresis on a 0.8% agarose gel.

| Product Test | ng/ μ l | 260/280 | 260/230 | Yield |
|-----------------|-------------|---------|---------|-------------|
| 1. Competitor Z | 162.5 | 2.00 | 2.07 | 8.1 μ g |
| 2. Competitor Z | 160.7 | 2.03 | 2.07 | 8.0 μ g |
| 3. Geneaid | 164.0 | 2.00 | 2.07 | 8.2 μ g |
| 4. Geneaid | 161.6 | 2.03 | 2.06 | 8.0 μ g |

GENEzol™ 96 Well TriRNA Pure Kit Real-Time PCR Data




Figure 2. Quantitative analysis of human beta globin mRNA extracted by GENEzol™ 96 Well TriRNA Pure Kit using a Taqman probe 1-step qRT-PCR assay. The assay was run on a BioRad IQ5 thermal cycler. The high yield, high quality extracted RNA was amplified quickly following a very short C_T (threshold cycle) compared to the RNA standards.

Related DNA/RNA Extraction Products

| RNA Purification | | |
|--|------------------|--------------------|
| Product | Package Size | Catalogue Number |
| Total RNA Mini Kit (Blood/Cultured Cell) | 50/100/300 preps | RB050/100/300 |
| Total RNA Mini Kit (Tissue) | 50/100/300 preps | RT050/100/300 |
| Total RNA Mini Kit (Plant) | 50/100/300 preps | RP050/100/300 |
| Presto™ Mini RNA Bacteria Kit | 50/100/300 preps | RBB050/100/300 |
| Presto™ Mini RNA Yeast Kit | 50/100/300 preps | RBY050/100/300 |
| miRNA Isolation Kit | 50/100 preps | RMI050/100 |
| GENEzol™ Reagent | 50/100/200 rxns | GZR050/100/200 |
| GENEzol™ TriRNA Bacteria Kit | 50/100 rxns | GZB050/100 |
| GENEzol™ TriRNA Pure Kit | 50/100/200 preps | GZX050/100/200 |
| GENEzol™ 96 Well TriRNA Pure Kit | 4/10 x 96 preps | 96GX04/10 |
| TriRNA Pure Kit | 50/100/200 preps | TRP050/100/200 |
| RNA Pure Kit | 50/100 preps | PR050/100 |
| Virus DNA/RNA Purification | | |
| Product | Package Size | Catalogue Number |
| Plant Virus RNA Kit | 50/100 preps | PVR050/100 |
| Viral Nucleic Acid Extraction Kit II | 50/100/300 preps | VR050/100/300 |
| Viral Nucleic Acid Extraction Kit III | 50/100/300 preps | VI050/100/300 |
| Genomic DNA Extraction | | |
| Product | Package Size | Catalogue Number |
| Genomic DNA Mini Kit (Blood/Cultured Cell) | 100/300 preps | GB100/300 |
| Genomic DNA Maxi Kit (Blood/Cultured Cell) | 10/25 preps | GDM010/25 |
| Genomic DNA Mini Kit (Tissue) | 50/100/300 preps | GT050/100/300 |
| gSYNC™ DNA Extraction Kit | 50/100/300 preps | GS050/100/300 |
| Genomic DNA Mini Kit (Plant) | 100 preps | GP100 |
| Genomic DNA Maxi Kit (Plant) | 10/25 preps | GPM010/25 |
| GENEzol™ DNA Reagent Plant | 100/200 rxns | GR100/200 |
| Presto™ Mini gDNA Yeast Kit | 100/300 preps | GBY100/300 |
| Presto™ Mini gDNA Bacteria Kit | 100/300 preps | GBB100/101/300/301 |
| Geneius™ Micro DNA Extraction Kit | 100/300 preps | GMB100/300 |
| Presto™ Buccal Swab gDNA Extraction Kit | 100/300 preps | GSK100/300 |
| Presto™ 96 Well Blood Genomic DNA Extraction Kit | 4/10 x 96 preps | 96GBP04/10 |
| Presto™ 96 Well Plant Genomic DNA Extraction Kit | 4/10 x 96 preps | 96GPP04/10 |

For additional product information please visit www.geneaid.com. Thank you!



Geneaid

www.geneaid.com



CERTIFICATE NO. QAIC/TW/50077