GENEzoI™ TriRNA Bacteria Kit

For research use only

Catalogue Numbers

GZB050, GZBD050 GZB100, GZBD100 GZB200, GZBD200

Quantity 50 rxns 100 rxns 200 rxns



Introduction

The GENEzol[™] TriRNA Bacteria Kit is a phenol and guanidine isothiocyanate plus spin column system for convenient purification of high-quality total RNA from bacteria samples. Bacterial cell walls are initially lysed using Lysozyme. The sample is then homogenized in GENEzol[™] Reagent without chloroform phase separation or isopropanol RNA precipitation. Following sample homogenization, simply bind, wash and elute the high-quality, total RNA in RNase-free Water and use in a variety of sensitive downstream applications.

Quality Control

The GENEzol™ TriRNA Bacteria Kit is tested on a lot-to-lot basis according to Geneaid's ISO-certified quality management system. 5 µl from a 50 µl eluate of purified RNA is analyzed by electrophoresis on a 1% agarose gel.

Advantages

- Purify total RNA within 20 minutes without chloroform phase separation or isopropanol RNA precipitation
- Up to: 1 x 10⁹ bacteria cells
- High quality RNA: A260/A280 >1.8, A260/A230 >1.8
- Applications: cDNA Library Construction, Cloning, RT-PCR (Endpoint), Real-Time PCR, Nuclease Protection Assays

Caution

GENEzol[™] Reagent contains phenol and guanidine isothiocyanate. During operation, always work in a fume hood, always wear a lab coat, disposable gloves, protective goggles and (anti-fog) procedure mask. Disposable/non-disposable glassware, plasticware and automatic pipettes should be sterile (RNase-free) and used only for RNA procedures.

Components and Storage

| Item | Volume | Product | Shipping | Storage |
|-------------------------------|---------------------------------|-------------|----------------------------|--------------------------------------|
| GENEzol™ Reagent | 4 ml | GZB004/D004 | | |
| | 40 ml | GZB050/D050 | room tomporature | |
| | 80 ml | GZB100/D100 | room temperature | dry at 2°C to 25°C |
| | 160 ml | GZB200/D200 | | |
| | 1.4 ml (0.6 ml) | GZB004/D004 | | dry at room temperature (15-25°C) |
| Pre-Wash Buffer ¹ | 21 ml (9 ml) | GZB050/D050 | room temperature | |
| (Add Ethanol) | 35 ml (15 ml) | GZB100/D100 | room temperature | |
| | 70 ml (30 ml) | GZB200/D200 | | |
| | 8 mg | GZB004/D004 | | -20°C |
| Lysozyme ² | 55 mg | GZB050/D050 | room temperature | |
| Lysozyme | 110 mg | GZB100/D100 | Toom temperature | -2010 |
| | 250 mg | GZB200/D200 | | |
| | 1.5 ml | GZB004/D004 | | |
| Bacteria Lysis Buffer | 15 ml | GZB050/D050 | room temperature | dry at room temperature |
| Bacteria Lysis Buller | 15 ml | GZB100/D100 | room temperature | (15-25°C) |
| | 30 ml | GZB200/D200 | 1 | |
| | 20 µl | GZBD004 | | |
| DNase I ³ (2U/µI) | 275 µl | GZBD050 |] | -20°C |
| Divase I [*] (20/µI) | 550 µl | GZBD100 | room temperature | -20°C |
| | 550 µl x 2 | GZBD200 | | |
| | 200 µl | GZBD004 | | |
| DNase I Reaction Buffer | 2.5 ml | GZBD050 | room tomporature | dry at room temperature (15-25°C) |
| Divase i Reaction Buller | 5 ml | GZBD100 | room temperature | |
| | 5 ml x 2 | GZBD200 | 1 | |
| | 2 ml (8 ml) | GZB004/D004 | | dry at room temperature |
| Wash Buffer ⁴ | 25 ml (100 ml) | GZB050/D050 | room temperature | |
| (Add Ethanol) | 50 ml (200 ml) | GZB100/D100 | room temperature | (15-25°C) |
| | 25 ml + 50 ml (100 ml + 200 ml) | GZB200/D200 | | |
| | 1 ml | GZB004/D004 | | dry at room temperature (15-25°C) |
| | 6 ml | GZB050/D050 | 1 | |
| RNase-free Water | 6 ml | GZB100/D100 | room temperature | |
| | 15 ml | GZB200/D200 | 1 | |
| | 4 | GZB004/D004 | | |
| RB Columns | 50 | GZB050/D050 | 1 | dry at room temperature |
| | 100 | GZB100/D100 | room temperature | (15-25°C) |
| | 200 | GZB200/D200 | 1 | · · · / |
| | 8 | GZB004/D004 | | |
| | 100 | GZB050/D050 | | dry at room temperature |
| 2 ml Collection Tubes | 200 | GZB100/D100 | room temperature (15-25°C) | |
| | 400 | GZB200/D200 | | () |

^{1,4}Add absolute ethanol (see the bottle label for volume) to Pre-Wash Buffer and Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

^{2, 3}DNase I and Lysozyme are shipped at room temperature and should be stored at -20°C for extended periods after receiving the kit.

RNA Purification Protocol Procedure

Please read the entire instruction manual prior to starting the Protocol Procedure.

Additional Requirements

absolute ethanol, 1.5 ml microcentrifuge tubes (RNase-free)

Optional Requirements

1 µL of 20 mM EGTA (pH=8.0) for Optional Step 2: DNA Digestion in Solution

1. Sample Homogenization and Lysis

Sample preparation should be performed at room temperature. To avoid DNA contamination of extracted RNA, be sure and use the indicated volume of GENEzol™ Reagent.

- 1. Transfer bacteria cells (up to 1 x 10⁹) to a 1.5 ml microcentrifuge tube (RNase-free).
- 2. Centrifuge at 12-16,000 x g for 2 minutes then remove the supernatant completely.
- 3. Weigh and transfer Lysozyme powder (1 mg/sample) to a new 1.5 ml microcentrifuge tube (RNase-free).
- 4. Add Bacteria Lysis Buffer (100 µl/sample) to the microcentrifuge tube containing Lysozyme.
- 5. Vortex the tube until the Lysozyme powder is completely dissolved.
- 6. Add 100 µl of Bacteria Lysis Buffer containing Lysozyme to the bacteria cell pellet.
- 7. Resuspend the cell pellet by vortex or pipetting.
- 8. Incubate the sample for 5 minutes at room temperature.
- 9. Add **700 µl of GENEzol™ Reagent**, mix well by pipette then incubate at room temperature for 5 minutes.

2. RNA Binding

- 1. Add 700 µl of absolute ethanol directly to the sample mixture.
- 2. Mix well by vortex then place a RB Column in a 2 ml Collection Tube.
- 3. Transfer 700 µl of the sample mixture to the RB Column.
- 4. Centrifuge at 14-16,000 x g for 1 minute then discard the flow-through.
- 5. Repeat the RNA Binding Step by transferring the remaining sample mixture to the RB Column.
- 6. Centrifuge at 14-16,000 x g for 1 minute then discard the flow-through.
- 7. Place the RB Column in a new 2 ml Collection Tube.

Optional Step 1: In Column DNase I Digestion IMPORTANT

DNA contamination is significantly reduced following In Column DNase I Digestion. However, traces of residual DNA may be detected in very sensitive applications. In this situation, please perform Optional Step 2: DNA Digestion In Solution instead to efficiently remove trace amounts of DNA. Standard DNase buffers are incompatible with In Column DNase I Digestion and may affect RNA integrity and reduce yield.

1. Add 400 µl of Wash Buffer (make sure ethanol was added) to the RB Column then centrifuge at 14-16,000 x g for 30 seconds.

- 2. Discard the flow-through and place the RB Column back in the 2 ml Collection Tube.
- 3. Prepare DNase I solution in a 1.5 ml microcentrifuge tube (RNase-free) as follows:

| DNase I | 5 μl (2 U/μl) |
|-------------------------|---------------|
| DNase I Reaction Buffer | 45 μl |
| Total volume | 50 μl |

4. Gently pipette the DNase I solution to mix (DO NOT vortex) then add DNase I solution (50 µl) into the CENTER of the RB column matrix.

5. Incubate the column for 15 minutes at room temperature (20-30°C) then proceed with RNA Wash.

3. RNA Wash

1. Add **400 µl of Pre-Wash Buffer (make sure ethanol was added)** to the **RB Column** then centrifuge at 14-16,000 x g for 30 seconds.

2. Discard the flow-through then place the **RB Column** back in the **2 ml Collection Tube**.

3. Add 600 µl of Wash Buffer (make sure ethanol was added) to the RB Column.

- 4. Centrifuge at 14-16,000 x g for 30 seconds. Discard the flow-through. Place the RB Column back in the 2 ml Collection Tube.
- 5. Add 600 µl of Wash Buffer (make sure ethanol was added) to the RB Column.
- 6. Centrifuge at 14-16,000 x g for 30 seconds then discard the flow-through.
- 7. Place the RB Column back in the 2 ml Collection Tube then centrifuge at 14-16,000 x g for 3 minutes to dry the column matrix.

4. RNA Elution

- 1. Place the dry RB Column in a clean 1.5 ml microcentrifuge tube (RNase-free).
- 2. Add 25-50 µl of RNase-free Water into the CENTER of the column matrix.
- 3. Let stand for at least 3 minutes to ensure the RNase-free Water is completely absorbed by the matrix.
- 4. Centrifuge at 14-16,000 x g for 1 minute to elute the purified RNA.

Optional Step 2: DNA Digestion In Solution

1. Prepare DNase I reaction in a 1.5 ml microcentrifuge tube (RNase-free) as follows:

| RNA in RNase-free water | 1-40 µl |
|-------------------------|----------------------------------|
| DNase I | 0.5 μl/μg RNA |
| DNase I Reaction Buffer | 5 μl |
| RNase-free water | add to final volume = 50 μ l |
| Total volume | 50 μΙ |

2. Gently pipette the DNase I reaction solution to mix (DO NOT vortex) then incubate the microcentrifuge tube at 37°C for 15-30 minutes.

3. Stop the reaction by adding 1 µl of 20 mM EGTA (pH=8.0) then incubate the microcentrifuge tube at 65°C for 10 minutes. NOTE: DNase I Reaction Buffer may cause aberrant migration or smearing of RNA on gels. If analyzing RNA by gel electrophoresis, repurify the RNA sample by using the Geneaid[™] RNA Cleanup Kit instead of stopping the reaction with EGTA.

GENEzol™ TriRNA Bacteria Kit Functional Test Data




Figure 1. RNA was extracted using the GENEzol[™] TriRNA Bacteria Kit. An *Escherichia coli* (1×10⁹) culture (OD600=2, 1 ml) was harvested by centrifugation at 16,000 x g for 2 minutes, followed by RNA extraction. 5 µl from a 50 µl eluate of RNA was analyzed by electrophoresis on a 1% agarose gel. M = Geneaid 1 Kb DNA Ladder

| Test | RNA Concentration | 260/280 | 260/230 | Yield |
|------|--------------------------|---------|---------|----------|
| 1 | 502.6 µg/ml | 2.06 | 2.18 | 25.13 μg |
| 2 | 518.3 μg/ml | 2.07 | 2.21 | 25.92 μg |
| 3 | 506.0 μg/ml | 2.08 | 2.24 | 25.30 μg |

M 1 2 3

Troubleshooting

| Problem | Cause | Solution |
|---|---|--|
| Low Yield | was incomplete B. Incorrect RNA elution C. Precipitates may form during the | A. Starting material should be reduced and completely dissolved in GENEzol[™] Reagent. B. Make sure RNase-free Water is added to the center of the RB Column and is absorbed completely. C. Reduce the sample amount to half of the original. |
| Degraded RNA | | A. Process or freeze samples immediately after collection.B. Extracted RNA should be stored at -70°C. |
| Low RNA A260/A280 | was insufficient for proper sample homogenization | A. Volume of GENEzol[™] Reagent is sample dependent and should be added according to the sample homogenization specifications. B. Wash the RB Column with ethanol added Wash Buffer 3 times. |
| Eluted RNA does not perform well in downstream applications | A. Residual ethanol contamination | A. Following the wash step, dry the RB Column with additional centrifugation at 14-16,000 x g for 5 minutes or incubate at 60°C for 5 minutes. |

Related DNA/RNA Purification and Extraction Products

| Plasmid DNA Purification | | |
|--|--------------------------------------|---------------------------------|
| Product | Package Size | Catalogue Number |
| Presto™ Mini Plasmid Kit | 100/300 preps | PDH100/300 |
| Presto™ Midi Plasmid Kit | 25 preps | PIF025 |
| Presto™ Midi Plasmid Kit (Endotoxin Free) | 25 preps | PIFE25 |
| Large Plasmid DNA Extraction Kit | 100/300 preps | PDL100/300 |
| Midiprep Spin Column Plasmid Kit | 25 preps | PA025 |
| Geneaid™ Plasmid Midi Kit | 25 preps | PI025 PIE25 |
| Geneaid™ Plasmid Midi Kit (Endotoxin Free) Presto™ Plasmid DNA Concentration Kit | 25 preps 250/500/1000 preps | PC0250/500/1000 |
| Geneaid™ Plasmid Maxi Kit | 10/25 preps | PM010/25 |
| Geneaid [™] Plasmid Maxi Kit (Endotoxin Free) | 10/25 preps | PME10/25 |
| Presto [™] 96 Well Plasmid Kit | 4/10 x 96 preps | 96PDV04/10, 96PDC04/10 |
| Presto [™] Plasmid 96 Well Binding Plate | 10 plates | 96PBP01 |
| Presto [™] Plasmid 96 Well FIlter Plate | 10 plates | 96PFP01 |
| Post Reaction DNA Purification | | |
| Product | Package Size | Catalogue Number |
| GenepHlow™ Gel Extraction Kit | 100/300 preps | DFG100/300 |
| GenepHlow [™] PCR Cleanup Kit | 100/300 preps | DFC100/300 |
| GenepHlow ™ Gel/PCR Kit | 100/300 preps | DFH100/300 |
| GenepHlow ™ DNA Cleanup Maxi Kit | 10/25 preps | DFM010/025 |
| Small DNA Fragments Extraction Kit | 100/300 preps | DF101/301 |
| Presto [™] Max Gel/PCR Kit (Large DNA Fragments) Presto [™] 96 Well PCR Cleanup Kit | 100/300 preps 4/10 x 96 preps | DFL100/300 96DFH04/10 |
| Presto™ 96 Well Gel Extraction Kit | 4/10 x 96 preps | 96DFR04/10 96DFG04/10 |
| Presto [™] PCR Cleanup Kit 96 Well Binding Plate | 10 plates | 96DBP01 |
| DNA Cleanup Kit | 100/300 preps | DP100/300 |
| G-25 Gel Filtration Desalting Column | 50 rxns | CG025 |
| G-50 Gel Filtration Dye Terminator Removal Column | 50 rxns | CG050 |
| 96-Well G-50 Gel Filtration Plate | 4/10 x 96 rxns | CGP04/10 |
| Gel Extraction Tool | 25 pcs | GXT025 |
| Genomic DNA Purification | | |
| Product | Package Size | Catalogue Number |
| Genomic DNA Mini Kit (Blood/Cultured Cell) | 100/300 preps | GB100/300 |
| Genomic DNA Midi Kit (Blood/Cultured Cell) | 25 preps | GDI25 |
| Genomic DNA Maxi Kit (Blood/Cultured Cell) | 10/25 preps | GDM10/25 |
| Genomic DNA Mini Kit (Tissue) | 50/100/300 preps | GT050/100/300 |
| gSYNC™ DNA Extraction Kit | 100/300 preps | GS100/300 GP100 |
| Genomic DNA Mini Kit (Plant) Genomic DNA Maxi Kit (Plant) | 100 preps 10/25 preps | GP100 GPM10/25 |
| Geneaid [™] DNA Isolation Kit (Blood) | 100/1,000 rxns | GEB100/01K(+) |
| Geneaid [™] DNA Isolation Kit (Bacteria) | 150/1,500 rxns | GEE150/1.5K(+) |
| Geneaid [™] DNA Isolation Kit (Tissue) | 150/1,500 rxns | GET150/1.5K(+) |
| Geneaid [™] DNA Isolation Kit (Cultured Cell) | 150/1,500 rxns | GEC150/1.5K(+) |
| GENEzol™ DNA Reagent Plant | 100/200 rxns | GR100/200 |
| Presto™ Mini gDNA Yeast Kit | 100/300 preps | GBY100/300 |
| Presto™ Mini gDNA Bacteria Kit | 100/300 preps | GBB100/101/300/301 |
| Geneius™ Micro DNA Extraction Kit Presto™ Buccal Swab gDNA Extraction Kit | 100/300 preps 100/300 preps | GMB100/300 GSK100/300 |
| Presto™ 96 Well Blood Genomic DNA Extraction Kit | 4/10 x 96 preps | 96GBP04/10 |
| DNA RNA Purification | 4/10 × 90 preps | 900BF04/10 |
| | Deelvere Circ | Catalagua Number |
| Product | Package Size | Catalogue Number |
| Presto [™] DNA/RNA/Protein Extraction Kit | 50/100 preps | DRP050/100 |
| Total RNA Purification | | |
| Product | Package Size | Catalogue Number |
| Total RNA Mini Kit (Blood/Cultured Cell) | 50/100/300 preps | RB050/100/300 |
| Total RNA Mini Kit (Tissue) | 50/100/300 preps | RT050/100/300 |
| Total RNA Mini Kit (Plant) Presto™ Mini RNA Bacteria Kit | 50/100/300 preps 50/100/300 preps | RP050/100/300 RBB050/100/300 |
| Presto™ Mini RNA Bacteria Kit | 50/100/300 preps | RBY050/100/300 |
| miRNA Isolation Kit | 50/100 preps | RMI050/100 |
| GENEzol™ Reagent | 100/200 rxns | GZR050/100/200 |
| GENEzol™ TriRNA Bacteria Kit | 50/100 rxns | GZB050/100 |
| GENEzol™ TriRNA Pure Kit | 50/100/200 preps | GZX050/100/200, GZXD050/100/200 |
| TriRNA Pure Kit | 50/100/200 preps | TRP050/100/200, TRPD050/100/200 |
| RNA Cleanup Kit | 50/100 preps | PR050/100 |
| GENEzol™ 96 Well TriRNA Pure Kit | 4/10 x 96 preps | 96GZX04/10 |